labCera®

DENTAL ZIRCONIA BLANK

KULLANIM KLAVUZU	TR
INSTRUCTIONS FOR USE	EN
GEBRAUCHSHINWEISE	DE
INSTRUCTIONS D'UTILISATION	FR
ISTRUZIONI PER L'USO	IT
INSTRUCCIONES DE USO	ES

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

EN Instructions for use

Safety notes

Please read the instructions for use carefully before taking the blank (made of zirconium oxide) out of its packaging. They contain important information which is essential to ensure both the correct processing and the safety of patients and users. Blanks made of Labcera Zrco are manufactured and tested according to highest quality standards.

Warning

When using this product, the dentist must consider possible interaction between this medical product and other products already in the patient's mouth.

General handling notes

Blanks made of Labcera Zrco are supplied in a pre-sintered White-fired condition. They have limited stability and strength and a residual porosity. Therefore, these blanks must be handled with care. The blanks must be stored in their original packaging at room temperature. They must not be subjected to impact or vibration. Contamination must absolutely be avoided. Make sure restorations are only handled with dry and clean hands or gloves, and that they are under no circumstances contaminated with liquids (such as adhesives or marker pens). Coolants will reduce the material's translucency.

Indication

Zirconia blanks used for manufacturing fixed dental prostheses. The material is suitable for single crowns and/or bridges consisting of up to 16* units. In the posterior region, the span between the abutments must not exceed two units. A cantilever bridge of premolar size is permissible.

Contraindications

In the event that there is an insufficient occlusal clearance and/or vertical prep wall, making the preparation unsuitable for an all-ceramic restoration, an alternative material must be chosen. Inlay bridges, endosseous implants and root posts are other contraindications. For bruxism, the vertical dimension must be observed. In these cases, please consult with the dentist.

Notes on construction

Wall thickness of the crowns in the sintered state must not be less than 0.5 mm. Connector geometry is of the greatest importance for crack resistance of abutment copings made of zirconium dioxide. Therefore, the cross section of the connectors must be as big as possible and should be a minimum of 9 mm² or, for cantilever bridges, 12 mm². For static reasons, the height of the connector is especially important.

Veneering:

According to manufacturer specifications, the veneering ceramic should be layered at a thickness between 0.7 and 1.5 mm, but < 2 mm. Furthermore, the abutment copings must be designed according to the general guidelines of digital dental technology.

Colouring the blanks

The restorations can be coloured with all approved coloring liquids according to manufacturer specifications labcera zirconia coloring liquids. Before the final sintering, the restorations which have been treated with liquids must always be dried under red light or in the furnace according to manufacturer specifications.

Processing

Only those machines and tools which are approved for processing pre-sintered blanks made of zirconium dioxide may be used for manufacturing restorations made of labcera zirconia blanks.


Caution! For multishade blanks the printed side is the occlusal side. Neither coolant nor compressed air should be used during processing. Please observe the instructions for use of your milling machine as well as the parameters of the CAD/CAM software.

Please separate the finished milled parts from the blanks carefully using suitable tools. Then, thin the edges which have been thickened by the machine and remove the tapping. Then, carry out a cut-back if this has not already been done in the CAD.

Final sintering

The final sintering, also called dense sintering, ensures that the milled restoration given its final properties. During the final sintering, the restoration shrinks by a precisely defined factor. In order to do this, the described procedure must be absolutely adhered to.

Sintering can take place in all common dental sintering furnaces which are approved for the sintering of restorations made of zirconium dioxide. Since zirconium dioxide is known to be a poor thermal conductor, it is recommended that you slowly heat the mounts to the required temperature (see sintering graph) and cool them down equally slowly. Sintering supports should be thinned out or milled as rods.

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Labcera Seramik San. Tic. Ltd. Sti